

TITLE: Evaluation of CRISPR-mediated nonhomologous end-joining (NHEJ) recombination

SOP Number: D-DEZ-PRO-015

Revision Number: 0

Effective Date: 07 Apr 2015

culture plate

2. The next day dilute 500ng of each plasmid into 200ul serum-free media per condition. Include a mock transfection control (200ul serum-free media only)
3. Dilute 5ul x (n+1) Lipofectamine 2000 into 100ul x (n+1) serum-free media.
4. Add 210ul Lipofectamine dilution to each diluted plasmid reaction (+ mock) mix.
5. Incubate at RT for at least 5 min.
6. While DNA and Lipofectamine are complexing, remove media from cells and replace with 500 ul complete media.
7. Add 200ul DNA/lipofectamine to appropriate wells and incubate overnight at 37°C, 5% CO₂.
8. The next day, remove complexes and replace with 1 ml complete media.

G. gDNA ISOLATION

1. 2-3 days post-transfection, isolate gDNA from each well.
2. Remove culture medium and rinse with 1X PBS.
3. Scrape cells into 1 ml of 1X PBS and centrifuge at 1000 rpm for 3 min.
4. Remove supernatant and isolate gDNA according to kit instructions – elute in 50ul sterile dH₂O.
5. Determine concentrations via Nanodrop and adjust to 100ng/ul.

H. SURVEYOR PCR

1. Set up PCR reactions as follows:

AccuPrime Reaction Buffer II 5ul
gDNA (100ng) 1ul
target oligo F (20nM) 0.5ul
target oligo R (20nM) 0.5ul
AccuPrime Taq 0.25ul
H₂O 42.75

2. Cycling parameters:

1x 94°C 3m
35x 94°C 15s
58°C 20s
72°C 1m
1x 72°C 5m

3. Check amplification by running 5ul on a 2% Egel. Estimate concentration of each sample by

TITLE: Evaluation of CRISPR-mediated nonhomologous end-joining (NHEJ) recombination

SOP Number: D-DEZ-PRO-015

Revision Number: 0

Effective Date: 07 Apr 2015

comparing to the 1000 bp band of the 100 bp NEB QuickLoad DNA ladder (95 ng). Amplicons should be at $\geq 20\text{ng}/\mu\text{l}$.

I. HETERODUPLEX FORMATION

1. Add 20ul each amplicon to new PCR tube.
2. Heteroduplex parameters:

95°C 10m
95-85°C -2°C/s
85°C 1m
85-75°C -0.3°C/s
75°C 1m
75-65°C -0.3°C/s
65°C 1m
65-55°C -0.3°C/s
55°C 1m
55-45°C -0.3°C/s
45°C 1m
45-35°C -0.3°C/s
35°C 1m
35-25°C -0.3°C/s
25°C 1m
25-4°C -0.3°C/s
4°C hold

J. SURVEYOR NUCLEASE DIGESTION

1. Set up digestions as follows:
 - heteroduplex DNA 20ul
 - MgCl₂ (150mM) 2.5ul
 - Surveyor S Enhancer 1ul
 - Surveyor S Nuclease 1ul

2. Incubate reaction for 1 hr at 42°C
3. Visualize digestions on 2% Egel.

K. QUANTIFICATION OF NHEJ

1. Dilute Surveyor PCR reactions 1:10 and subclone 4ul into pCR2.1-TOPO-TA vector.
2. Transform 4ul TOPO TA reactions into TOP10 OneShots.
3. Plate 150ul outgrowth onto LB + AIX plates.
4. Submit plates to Functional Biosciences for prep and sequencing with M13 (-20) Forward primer.

TITLE: Evaluation of CRISPR-mediated nonhomologous end-joining (NHEJ) recombination

SOP Number: D-DEZ-PRO-015

Revision Number: 0

Effective Date: 07 Apr 2015

L. * TARGET OLIGO DESIGN

1. When designing target oligos use PrimerQuest on IDT website
<http://www.idtdna.com/Primerquest/Home/Index>
2. Choose ~400 bp on either side of gRNA target sequences to input.
3. Choose primer pairs with amplicons of ~400 bp to ~600 bp.
4. Avoid pairs with gRNA target sequences at center.

M. HISTORY

Effective Date	Revision	Change
07 Apr 2015	0	Initial issue of SOP